


Satelec News • January 2010 • n°1



New Tip's blister !

Since January 1st 2010, Satelec has developed new blister design to highlight its original tips.

Front's Blister



Sticker for Tip number and Color Coding System

Guarantee of Satelec Original Tip

Back's Blister



Cleaning and sterilization Procedure

Tracability number



Pull the tag


Easy Opening on back's blister

See overside for packaging recommendations.





The tips under new blisters have to be packaged with cautions when being store or being send to your customers
Please keep on mind this packaging's recommendations.



Put tips 2 per 2, back to back

Shipment for 3 tips + wrench



Shipment for 6 tips
And Shipment for 10 tips



Box Reference :

- 3 tips + wrench / 6 tips : N12007
- 10 tips : N08060

